MANUAL | MANUEL | MANUALE

hyghtec® basic-plus

Gebrauchsanweisung | Instructions for use
Mode d'emploi | Istruzioni per l'uso

Deutsch | English | Français | Italiano

Deutsch

hyghtec®

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English

hyghtec®

... innovative catheter design with polyurethane

hyghtec stool drainage systems are based on the innovative combination of micro-thin-walled balloon films & flexible tube components made from polyurethane (PUR).

hyghtec stool drainage systems optimally adapt to fit the individual anatomy and function.

hyghtec stool drainage systems absorb the forces at work in the rectum and use them smartly for efficient trans-anal sealing.

hyghte products offer solutions for all indications:

› **hyghte basic-plus**

... defines effective bowel management
in normal and intensive care settings

› **hyghte basic-ped**

... facilitates closed-system stool drainage
even for paediatric patients

› **hyghte basic-home**

... supports home care routines

› **hyghte i+d**

... for intensive, colorectal lavage with
subsequent stool drainage

› **hyghte irrigator**


... for leak-proof, therapeutic and diagnostic
irrigation

hygtec® basic-plus

DESIGN AND FUNCTION

... a functionally segmented balloon attached to a flexible, collapsible tube

- (A) Balloon segment that anchors the system intra-rectally
- (B) Leak-proof, trans-anal balloon segment
- (C) Shaft tube
- (D) Gel tip
- (E) Coloured ring for checking correct positioning
- (F) Stool irrigation tube

Leak-proof balloon segment for trans-anal anchorage

-) Micro-thin & tear-resistant
-) Fully moulded
-) Inserted semi-deflated at low inflation pressure


Shaft tube

-) Flexible & self-expanding
-) Can be folded in the axial and radial directions without being damaged

hyghtec absorbs the forces at work in the rectum via the intra-rectal segment (A) and uses them in the trans-anal segment (B) for sealing purposes .

(G) Irrigation channel

(H) Irrigation aperture directed into the cone

System irrigation

- › Cleaning of stool-draining surfaces through regular system irrigation
- › Removal of stool from drainage aperture by releasing the irrigation fluid into the interior of the tip

Rectal irrigation

- › Rectal enema can be performed by closing off the stool-draining tube using the easy-to-use clamping strap


hyghtec is fitted with an effective odour barrier. Odours are substantially reduced even when hyghtec is used over a prolonged period.

hygtec® basic-plus

QUICK GUIDE

INSERTION

- 1** Digital inspection of anus and rectum
- 2** Deflate balloon completely
- 3** Apply lubricant to the first third of the balloon
- 4** Insert the first two-thirds of the balloon into the anus
 -) Hold the tip like a suppository when inserting

- 5** Fill with **85 ml of air**
- 6** Check the trans-anal position
 -) Yellow marking outside of the anus

USE

7

Once per nursing shift:

- › Check position
- › Check balloon inflation
- › Check flow rate

8

Once a day: System irrigation

9

Optional: rectal irrigation

10

After 48 hours: change the stool collection bag
(recommendation)

11

After 72 hours: completely deflate and
re-inflate the balloon

12


Maximum duration of use: < 30 days

SYSTEM IRRIGATION

- 1** Check position and balloon inflation
- 2** Check the temperature of the irrigation fluid (body temperature)
- 3** Smooth out the irrigation tubing to the stool collection bag

SAMPLING


- 1** Clamp off the stool-draining tubing by tightening the clamping strap below the stool removal port
- 2** Insert the hub of the sampling syringe directly into the sampling port

- 3** After taking the sample, loosen the clamping strap fully and position it at the stool collection bag

RECTAL IRRIGATION

- 1 Check position and balloon inflation
- 2 Check the temperature of the irrigation fluid (body temperature)
- 3 Stop the drainage by tightening the clamping strap approximately 10 cm from the anus

- 4 After irrigation, loosen the clamping strap fully and position it at the stool collection bag

hyghtec® basic-plus

INSTRUCTIONS FOR USE

Product description:

hyghtec basic-plus stool drainage system

Product design and function:

Unlike conventional stool drainage systems, hyghtec basic-plus is fitted with a micro-thin, hourglass-shaped balloon that positions the drainage system past the sphincter muscle and seals it. Inside the balloon is a stool-draining shaft tube which is stabilised by its special ridged design. It has crumple-like zones which make it compressible. The movements of the shaft components directed towards the patient are therefore buffered. With a normal sphincter tone, the shaft tube in the anus collapses to a small, non-irritating diameter. When the sphincter relaxes, the radially folded shaft tube expands and opens. The trans-anally positioned balloon has already been moulded to its full working dimensions during production. It does not, therefore, have to be stretched under pressure. Unlike conventional silicone drainage systems, the balloon is inflated with air, not water. This quantity of air is used so that the balloon remains flaccid and is not fully inflated. The thin, flaccid, membranous balloon sheath optimally conforms to the surfaces of the rectum and anus. This ensures that efficient sealing is achieved even at the lowest inflation pressure and in the case of watery stools. The balloon is inflated with the hyghtec inflating syringe supplied. The syringe is graduated and has on it instructions on the exclusive use of air as the filling medium. Flexible irrigation tubing is connected to the trans-anally positioned top section of the hyghtec basic-plus. The bayonet lock at the lower end of the irrigation tubing facilitates replacement of the stool collection bag (hyghtec bag). The bag is fitted with an odour filter to vent the drainage system. The bag contents are concealed from the user and can be viewed through the transparent rear side of the stool collection bag. The hyghtec basic-plus device has an integrated supply line for rinsing and irrigation fluid. hyghtec stool drainage systems are fitted with an effective odour barrier. Odours are effectively contained even with prolonged use.

Product indication and contraindications:

Faecal excretion is often disrupted in critically ill patients. There are periods with constipation as well as periods with thin watery stools. hyghtec basic-plus is designed specifically for comprehensive bowel management. hyghtec's highly efficient trans-anal

seal ensures intensive, leak-free irrigation as well as a reliable seal for continuous, therapeutic stool liquefaction.

The hygtec basic-plus drainage system must be prescribed by a doctor.

The indication for the product as follows:

-) Continuous, trans-anal drainage of liquid stools for <30 days

Contraindications:

-) Formed stools that will not flow easily
-) Acute traumatic injury to the rectum or anus
-) Previous colorectal surgery within the last 12 months
-) Rectal or anal strictures, stenoses or tumours that pose an obstacle to the trans-anal use and positioning of the drainage system
-) A narrow rectum with insufficient space for the intrarectally anchored drainage system balloon
-) Inflammatory degenerative changes in the rectal mucosa such as severe proctitis or ischaemic proctitis
-) Advanced haemorrhoids prone to bleeding
-) Patients with spinal cord injuries or rectal, anal or pre-anal hyposensitivity
-) Spinal shock phase
-) Circulatory therapy requiring catecholamines
-) Children less than 12 years of age
-) Pregnancy
-) Known allergies to polyurethane or PVC
-) Use of the drainage system >30 days

All contraindications should be considered on a case-by-case basis by the supervising physician.

General risks associated with the continuous use of the stool drainage system:

-) Irritation due to pressure or friction, inflammation or ischaemic-necrotic damage to the rectum and/or anus
-) Transient or permanent functional damage to the anal sphincter muscle
-) Ischaemic and/or mechanical perforations of the intestinal wall
-) Rinsing or irrigation with media used at the wrong temperature, media that can cause damage to the mucosa or which otherwise cannot be tolerated

Specific risks associated with the use of this product:

-) Balloon inflation with a filling volume exceeding 85 ml of air or permanent over-inflation of the balloon
-) Use of a liquid medium to fill the balloon

Note the following warnings and precautions to avoid side effects and misuse:

-) Only use the product if the packaging is intact.
-) It is imperative to perform a digital inspection of the anus and rectum before inserting the hyghtec basic-plus.
-) The rectum must have enough room for the anterior, rectally positioned segment of the balloon to unfold. The rectum must be completely free from stools. If stool residues are still visible in the rectal ampulla, a cleansing enema should be performed before inserting the drainage system.
-) If the drainage system is in the correct trans-anal position, the yellow coloured ring will be visible outside the anus.
-) The balloon segment protruding from the anus can be palpated to check the inflation of the anchoring/sealing balloon. If the pre-anal balloon remains flaccid after being inflated with 85 ml of air, the drainage system must be removed and replaced.
-) Avoid constant tugging on the irrigation tubing in the axial direction.
-) The drainage system must not be used in conjunction with other catheters or devices positioned in the anus.
-) Examine the irrigated stools for any signs of blood or bloody secretions. If any bleeding is observed during use, the patient must be examined immediately for signs of ischaemic-necrotic tissue changes (rule out pressure damage).
-) Check for abdominal bloating when using the drainage system. Check intestinal activity using auscultation at least once per nursing shift.
-) Inform the treating physician immediately if the patient experiences pain or irritation in the rectum or anus.
-) Solid or formed stools will not pass through the drainage system. The use of hyghtec basic-plus should be discontinued if the stools are not a liquid enough consistency.**
-) Completely deflate the balloon before removing the product.
-) Ensure that only body-temperature fluids are inserted into the patient's rectum. If chemically or pharmacologically active additives are used, these should be initially checked for any interactions with the basic materials used in the drainage system (PUR inside the body, PVC outside the body). Always discuss the potential administration of active substances with the responsible physician or pharmacist. Clean the irrigation channel with 50 ml of saline solution before and after administering active substances.
-) Note any shift in the body's fluid retention during rinsing and irrigation, especially in patients with circulatory and renal insufficiency.
-) Use the device and accessories on no more than one individual patient.

Inserting and using the drainage system:

-) Connect the drainage system to the stool collection bag (hyghtec bag).
-) Attach the stool collection bag by its strap to the lower bed frame and uncoil the irrigation tube.

-) Deflate the balloon completely before inserting the drainage system.
-) Moisten the front third of the balloon with a lubricant (water-soluble substances are recommended).
-) Hold the green insertion tip with the thumb, index and middle finger, and insert approximately 2/3 of the length of the balloon into the patient's anus.
-) **Inflate the balloon with 85 ml of air using the filling syringe supplied. Do not exceed the filling quantity. Only use air.**
-) Position the stool-draining tubing as close to the patient as possible to avoid any twisting or kinking of the irrigation tubing. Ensure that the patient is not lying on the irrigation tubing.
-) Check the trans-anal position of the drainage system at least once during each nursing shift. If the device has been positioned correctly, the yellow coloured ring will be fully visible just outside of the anus.
-) If the drainage system is expelled, the balloon must be deflated completely before being reinserted. Once inserted, the balloon is inflated with a single syringe filling of 85 ml of air.
-) Check that the stools are unobstructed and free-flowing at least **once per nursing shift**.
-) Clean the drainage system at least **once a day** by irrigating the system.
-) Check the balloon inflation status at least once a day. To do this, inspect the balloon segment protruding from the anus.
-) Completely deflate the balloon **every 72 hours** and re-inflate with 85 ml of air.
-) For rectal irrigation, the tubing is clamped off for the duration of the enema using the clamping strap, located approximately a hand's width from the anus. Once opened, the strap is positioned near the stool collection bag.
-) hyghtec basic-plus is fitted with a port for the removal of stool samples. The sampling port is located in the lower third of the stool-draining tubing. To take a stool sample, the stool can be contained by clamping the clamping strap below the port.
-) Regardless of how full it is, change the stool collection bag **after 48 hours** (recommendation).
-) The hyghtec basic-plus drainage system must not be used for more than 30 days.

Product disposal:

The used, contaminated system must be hygienically disposed of on removal from the rectum. hyghtec drainage systems must never be cleaned and reused.

hyghtec® basic-plus

BLISTER PACK

... designed for normal and intensive care settings

Contents of the blister pack:

- › 1 stool drainage system
- › 1 stool collection bag (hyghtec bag)
- › Graduated inflation syringe
- › Clamping strap
- › Manual (Instructions for Use)

V01-10017 – selling unit: 1 drainage system

hyghte® bag

... the stool collection bag for all hyghte stool drainage systems

Features:

- › 2 litre capacity
- › Odour-absorbing ventilating filter
- › Concealed front view
- › Transparent rear view
- › Leak-proof bayonet lock

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Creative Balloons GmbH