

AnapnoGuard™

Innovative Airway Management Solution



HOSPITECH
RESPIRATION

The AnapnoGuard CUFFILL

Accurate Cuff Pressure Control Made Easy

The compact and intuitive AG CUFFILL device is the most accurate solution for measuring both pressure and volume of airway cuffs in all clinical settings. Its pocket-size, syringe-like design allows a simple and easy operation by medical professionals, including first-responders and hospital staff, and by users in the home care environment.

The AG CUFFILL is designed to enhance patient safety by minimizing the potential of human error while measuring and inflating airway cuffs, as well as diminishing the risk of cross-contamination, when used for a single patient.

It is recommended for single-patient use in order to reduce the risk of cross-contamination.



AnapnoGuard™

Innovative Airway Management Solution



HOSPITECH
RESPIRATION

AG CUFFILL Highlights

- Volume 10 cc
- Measurement range 0-99 (mmHg/cmH₂O)
- Measurement accuracy ± 2 (mmHg/cmH₂O)
- Disposable, enables 100 measurements
- FDA and CE approved

Scale	Cat No.
-------	---------

mmHg	HSCUFF0031
------	------------

cmH ₂ O	HSCUFF0041
--------------------	------------

Instructions for Use

- Turn the CUFFILL device ON by pressing the yellow button (it will stay ON for 60 seconds)

For cuff pressure measurement:

- The CUFFILL device should be connected to the tube when the plunger is pushed to the distal end of the syringe throughout the entire measurement

For cuff inflation/deflation

- The CUFFILL device should be connected to the tube when the plunger is pulled back about 3 cc (or as required)
- The plunger should be pushed/pulled until the required cuff pressure is presented



- Cuff pressure will change with the breathing cycle
- Cuff pressure to be recorded is the lowest pressure presented

About Hospitech

Hospitech Respiration develops innovative airway management solutions to improve patient safety and reduce complications and hospitalization costs of mechanically ventilated patients